



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Auto Refractometer

Item Code:

#### Product Description

- 1, Faster and more accurate measurement: measurement time is reduced half compare with a conventional product
- 2, More elegant and tailor-made design .
- 3, Measurement minimum pupil diameter -2.0mm: pupils of patient all over the world can be measured easily
- 4, Selectable functions mode: you can select and change diopter steps, vertex distance, cylinder mode, measuring mode, shop name ,date and time, The data measured can be used for different modes of STEP, VD and CYL

#### Product Feature

- 1, Measuring Range:  
Sphere -20~+20D (VD=12) 0.125D steps  
Cylinder -8~+8D 0.12D/0.125D steps  
Axis 0~180°?1°steps
- 2, Pupil Distance 45~85mm, 1mm steps
- 3, Vertex Distance 0mm, 12mm, 13.75mm
- 4, Min. Pupil Size  $\varnothing$ 2.0mm
- 5, Chart Auto fog
- 6, Display 5" CRT Monitor
- 7, Print Thermal Printer
- 8, Power Supply 100~120v/200~240V?  
50/60HZ.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : rajasmedical@gmail.com  
www.rajasingia.com