



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Autoclave Horizontal Rectangular

Item Code: S.T.R-01

our Autoclave Horizontal Rectangular High Pressure. This is a triple walled autoclave used to sterilize the medical and testing equipment at specific temperature. Our high pressure steam sterilizers are highly effective in cleaning hospital appliances from the infectious organisms and chemicals. The entire array is made of premium grades raw materials supplied to us by the prominent vendors of the industry. We supply our comprehensive array at affordable rates.

Technical Specifications:

Construction	Sterilizer has inner and outer wall made of 304 grade stainless steel, middle jacket made of stainless steel or brass and gun metal.
Technical Data	Steam Working Pressure 1.26 kgf/cm ² (2.2 kgf/cm ² in case High speed sterilizer
Working Temperature	121oC
Operating Voltage	240 volts single phase, AC supply, 50 Hz.
Sterilization Period	45 to 50 minutes
Components	Radiant locking system, cold water drainage system
Steam Exhaust	5 to 7 provided with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the pressure inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber made of stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level indicator to see the water inside, heavy duty power plug with socket provided with the unit
Accessories	One spare element, and two gasket , our Regular sizes are 900X500, and 1050X750. Other sizes also available on request.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com