



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Autoclave Kerosene Type Sterilizer

Item Code: S.T.R-10

Rajas Autoclave Non Electric stove driven made of full stainless steel and with heavy gun metal lid , the portable autoclave used for the place where electric supply is not available, the equipment provide best sterilization at 120 Degree , one pressure gauge, one safety whistle , one safety valve, wing nuts to lock the lid, handle provided to carry the lid, best of it's class total portable available in differnt sizes as per customer requirments.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com