



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Autoclave Table Top Microprocessor Control

**Item Code:**

Fully automatic, instant Autoclave of 20 ltr capacity, having complete body of stainless steel of grade SS 304. It should provide sterilization at 121C and 134C with pre selectable sterilization time. It have pre vacuum cycle to remove air from the chamber. Water tank inbuilt with automatic water level sensors. It add water in the chamber automatically. It have automatic pressure door locking system when sterilization under pressure is in progress, to prevent accidental door opening. Automatic process termination in case, door is opened accidentally. It have dry cycle with selectable time period. Audio visual alarm to indicate end of process. Compatible sterilization tray be given with the autoclave. Visual indicator provided by Red LED to indicate process status. Safety valve to release excess pressure to prevent accident. Alpha numeric display to set the protocol with soft keypad.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)