



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Basin On Stand Double

Item Code: H.F-11

basin on stand Double. Frame made of steel construction duly powder coated paint finish, two rings fixed on the top of the steel rod duly supported with steel strip, two basin of stainless steel duly bended from corners, Base made of Poly urethane material with five castors Black in colour mobile unit.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com