



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Bed Side Screen Four Fold

Item Code:

Bed Side Screen with Four folds. Frame of Bed side Screen made from M.S. pipe which is duly powder coated paint finish Hospital Bed Side Screen, 4-Fold Size : 244cm L x 165 cm H (opened) : 60 cm L 165cm (closed) Hospital Bed side screen 4 Fold has Removable Feet which allow a compact packing for easy transportation. the Unit is on wheel for easy placement in the room or ward.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com