



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Blood Bag Tube Sealer

Item Code: BB-04

Rajas Blood Bag tube Sealer equipment to seal the blood bag tube by high frequency dielectric heating system without causing haemolysis, leakage of blood. The sealing can be achieved in less than 1.5 seconds per tube by simple operation. Complete sealing of the blood bag tube is annexure by high frequency dielectric heating. Sealing of each tube takes approx 1.5 seconds and is thus very efficient. Pulling the sealing section to both sides can readily cut off sealed tube Input Voltage 230 V AC 230 V AC(12 V SMPS) Radio Frequency 40.68 MHZ 40.68 MHZ Sealing time Approx 1.5 sec. Approx 1.5 sec. Power Consumption 100 V A 100 V A Weight 4.4 Kg approx 5.4 Kg approx Size 170 x 210 x 110mm (WxDxH) 175 x 250 x 110mm (WxDxH) Battery Backup 150 sealing.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)