



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Blood Bank Refrigerator

Item Code: BB-01

Our blood bank refrigerators are mounted on a sturdy angle frame the inner chamber is made of stainless steel and outer chamber made of M.S.duly powder coated paint the temperature is maintained with 4 degree to 6 degree for storing blood bottles or bags temp is controlled by Digital indicator cum controller and triggers the audio visual alarm in the event of fluctuation exceeding permitted tolerance unit is also fitted with door operated alarm. And Our blood bank refrigerators have power supply failure alarm as a safety device. Our blood bank refrigerators drawer are provided for holding the bags to work 220-230V and 50 to 60HZ the accuracy level of temperature is set even 0.5 degree or may be more precise, the blood bank refrigerated fitted with data logger which is connected with your PC and with RS 232 port to maintain all the data on your PC for temperature record, this latest technique help to maintain the date on computer and avoid hassle to change the paper and pen in thermograph and then to keep thermograph paper records, we also offer many customize feature with our blood bank on the demand of our customer with additional cost.( thermograph is also available on request)

#### Features:

- Temperature is maintained with 4 degree to 6 degree Celsius
- Digital indicator cum controller for controlling the temperature
- Audio visual alarm in the event of fluctuation exceeding
- Tolerance unit is fitted with door operated alarm
- Power supply failure alarm
- Data logging facility for 100 readings
- Alpha numeric LCD Display
- Capacity: 50, 110, 200, 250 and 360 blood bags

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : rajasmedical@gmail.com  
www.rajasindia.com