



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Blood Collection Monitor

Item Code: BB-03

Customers can avail from us a wide range of **blood collection monitors (item code: b. B-03)** that are designed using latest technology. This compact machine ensures the anticoagulant and mixing of blood without harming blood cells during the process of collecting blood from the donor. Our designed monitors are available with numerous innovative features like pausing during the process of blood collection, easily set the volume to 1 ml increment, preset volume can be changed, gram and totaling of time during the process, an audio video alarm that indicates flawless flow and rings if the flow decreases to 15 ml or low and can also rings after accomplishment of the collection process. Furthermore, a micro controller program is installed into our system to ensure high accuracy and to stop the excess blood flow from the required volume; an auto tare facility motorized pinching system is also available in our machine.

Features:

- Trouble free operation
- Notification of end of the process
- High accuracy

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com