



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Blood Donation Chair Cum Bed Portable

Item Code: BB-24

Rajas Blood Donation chair portable and folding type easy to carry for outdoor activity especially for Blood Donation camps and other outdoor activity. the chair can handle weight up to 220 kg construction of Aluminum pipe back of the chair is 50x70 cms and seat of the chair is 50x100 cms length of the chair is from 170 to 190 cms. the chair has mattress which can be removed with the thickness of 4-5 cms which is washable and easy to clean fabricated with semi leather material. the chair has provision to set two armrest which is removable and one head rest on head side. total weight of the chair is 10 to 15 kg.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com