



Blood Donor Bed (Stationary)

Item Code: BB-23

Donor bed for use in a stationary setting For donor rest area. Height approx 50-80 cm. Built for easy transportation, ie. Foldable & light weight construction. Collapsible. able to carry 150 kg weight. Length max 6 approx feet the Donour bed have both side arm rest which has height adjustment and can rotate at 360 degree as per need of the user, one side there is one tray for holding collection monitor or blood bag, head side of the bed is adjustable on desirable position with help of knob provided on the back side, uphostely duly made by semi leather type material which is easy to clean, frame of the bed made of square pipe of M.S. duly powder coated paint finish.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com