



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Blood Donor couch

Item Code: BB-20

Blood Donor Couch for mobile blood collection. For donor rest area. Height approx. 50 cm. Built for easy transportation, Foldable and lightweight construction of round pipe of Stainless Steel. Collapsible. able to carry 150 kg weight. Length max 5-6 feet this is chair cum couch easy to fold and easy to open as couch, foam padded seat, back and foot side, all the three parts are movable as per need of the user, stainless steel construction, both the arms are removable as per user need, arm are rotatable as well on 360 degree, foam padded armrest allow full comfort to the patient while donating the blood. The couch allow Trendelenburg position and reverse Trendelenburg position

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)