



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Blood Donor Couch (Mobile)

Item Code: BB-22

Couch, donor for mobile sessions. frame made of steel round pipe duly powder coated finish. For donor rest area. Height approx.50 cm. Built for easy transportation , ie. Foldable and light weight construction. Collapsible. Able to carry 150 kg weight. Length max 5-6 feet this is folding donor chair easy to fold and easy to open as needed, comfortable seat, back and foot side, all the three parts are movable as per need of the user, steel construction, both the arms are foam padded arm rest allow full comfort to the patient while donating the blood. Seat back and foot section is covered with tetron fabric which is easily washable and enough sturdy to hold the weight of the patient

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasingia.com](http://www.rajasingia.com)