



Blood Donor Couch (Stationary)

Item Code: BB-21

Stationary Couch cum donor chair for stationary use For Blood donor rest area. Height approx.50 cm. Built for easy transportation, ie. light weight construction. frame made of steel construction duly powder coated paint finish, may be folding. able to carry 150kg weight. Length max 5-6 feet this is chair very comfortable back side of the chair is adjustable with mechanish, both the arms are removable and rotatable as well , easy to clean seat and back, stationary chair can be transport from one to another place

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com