



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Cold Box

Item Code: BB-17

For transport of blood. Heavy duty, domestic type. Approx. 60-70 cm x 40-45 cm x 42-46cm. Able to contain 20 blood bags and appropriate cold packs that can fit into its compartment. Able to maintain +2 to +8 degree C for 12 hours. Have 2 handles to enable one or two persons to carry. Construction of PVC material easy to handle for outdoor activity, best A class insulataion for best result

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com