



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Crash Cart

Item Code:

Crash Cart made of Best quality heavy gauge tubular pipe and metal sheet construction, the unit moves on the four castors in which two are with break and two are without breaks, Top and lower shelves of the unit made of stainless steel, Six number of total drawers provided on the unit three on the each sides made of plastic and the plastic tray of different colors on the frame to keep the emergency medicines the top of the unit covered from all the three sides with the help of the stainless steel rod, the handle provided on the unit to move around. The lower shelves also covered from all the three sides with the help for the stainless steel rods for the safety of the material. I.V rod and cage for oxygen cylinder provided in the trolley.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com