



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Deep Freezers

Item Code: BB-10

Rajas micro controller Based Cryoprecipitate freezer, inner chamber made of heavy gauge of stainless steel (304) outer body of mild steel heavy gauge duty powder coated paint finish, temp up to -20C -40C and -80C( as per requirements) with built in Micro controlled based temperature indicator cum controller having audio visual alarm for lo and high temperature. Heavy Duty Stainless Steel Basket. Display of Set and Process value temprature Extra-strength non CFC refrigeration System. With Danfoss Compressor. unit is supplied with tow option one is with data logger which allow to keep all the data in your PC and supported with compatiabile software, another option is chart recorder. customer can choose any of the option as per choice The unit is mounted on Heavy duty Castor wheel (Lockable) for easy mobility. 120 mm thick CFC Free PUF (Poly Urethane Foam) Insulation to ensure temperature stability and save power consumption. Stainless Steel inner door with insulation. Double magnetic gasket in main door. Available in 80 & 135 Bags Capacity. Available in -20, -40 and -80 Degree.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)