



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Dental Chair Hydraulic (Manual)

Item Code:

Dental chair manually operated hydraulic chair with 'Z' base, chair has all positions such as up-down, back complete flat and up. Headrest rotates to 180 Degree and could used as pediatric head rest. Chair can be fitted manually operated unit or fully automatic unit as desired. Surgical stool, height adjusted by pneumatic piston with adjustable back rest. Provision of Inbuilt suction unit one three way syringe and provision for air motor and air rotor hand piece, Halogen light having two intensity of 9 V and 12 Volt, spittoon and tumbler operated by solenoid valve a x-ray view box provided with the unit to view the x ray on site, water attachment with water tank provided, heavy base sheet finish with off white powder coated paint and the color of the seat as per choice of the customer a thick lamination sheet provided on the seat to avoid the damage due to any liquid substance during procedure foot control hydraulic system given smooth movement to the chair for up and down.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com