



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Dental Compressor Oil Free

Item Code: DE-03

Rajas 1/2 HP Dental compressor seal unit fitted with safety valve for air release excess sorted, auto cut off system, pressure gauge and drain valve air rotor air compressor mounted on high pressure tank entirely powder coated paint, Power ½ HP Amperes 4.5 Flow 40 to 50 Liter per minute Speed 1000 rpm Noise 50 to 60 dB (A) Power supply 220V, 50 Hz Tank 25L Pressure 350 MMHG Weight 35 to 37kgs

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com