



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Dental X Ray Machine

**Item Code:**

The dental X-ray package is a true insurance of quality with all the benefits of the RAJAS technology. The technical features of our tube heads make them the safest on the market. An extremely high dielectric coefficient is obtained through use of specially treated oil. The high voltage transformer comes with the best insulators. special expansion chamber compensates for the dilation of the oil caused by heat from the tube. The protection against leakage radiation is stricter than the strictest international standards. The hooded tube with grid compensates for voltage fluctuation. Input voltage: 230V AC single phase Timer: Electric with digital display of mAs Output: 8 mA - 70 KVP Stand: Floor model/Wall model-ray tube: 1.5mm focal Spot Cone: Paralle

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)