



An ISO 9001 & 14001 Co.

# HOSPITAL & MEDICAL EQUIPMENT



## Dimroth Condenser for Soxhlet

Item Code:

**Dimroth Condenser for Soxhlet** as per DIN Standard made from borosilicate glass 3.3 (Supplied with either glass or polypropylene side arm which has screw thread connector for easy fitting on the tubing)

Cat. No.	Extractor Capacity ml	Cone
R55.030.29/32	30	29/32
R55.070.34/35	70	34/35
R55.100/250.45/	100/250	45/40

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)