




An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Scalp Vein Set

Item Code:

Disposable Scalp vein set presented in single air tight packing, sterilized by gamma radiation, available in 16 G to 24 G adequate long tube with locking at the end, butterfly wings in different colour, sharp edge needle made of high grade stainless steel.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)