



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT

Doctor Stool

Item Code: DE-02

With immense experience in the field, we offer a broad array of **Doctor Stool** that are highly appreciated in the medical industry. These stools are manufactured from top quality raw material using most advanced technologies. Our **pneumatic surgical stools** are designed as per the requirement of modern and hi-tech medical organizations. A pneumatic piston is provided to adjust the height as required with adjustable backrest to offer more ease to the dentist as well as the patient. Cushioned Back and seat with Back rest, Hydraulic up and Down Movement, five spokes base Chrome plated. Height adjustable with help of lever under the seat.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com