



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Dressing Trolley with SS Bowl & Bucket

**Item Code:**

Dressing trolley made of Stainless steel construction total overall size of 30"x20"x32" the trolley with S S. Top and frame also made of stainless steel pipe round type units moved on castors two with brakes and two without breaks provided with two S. S. Shelves, hinges of stainless steel rod on the all the three sides for the support on both the shelves upper and lower. Utensils provided made of stainless steel and steel rod ring to hold the same properly welded with the trolley.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)