



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



ECT Machines

Item Code:

Our ECT Machines delivers pulsed, bi-directional square wave stimulus. The current is constant at 0.9 ampere and the voltage varies automatically by a built in microcomputer as per the impedance of the patient. We conceived and designed the rajas keeping three major considerations in mind i. E. Its stimulus parameters had to be chosen on the basis of research data. (Maxwell 1968, weaver teal 1977 wiener 1980). It should have the capacity to produce seizure virtually 100%. It should be safe and easy to operate. Parameters display on LCD. Auto stimulus voltage. Auto impedance check. Output display in joules as well as in milli coulombs. In built rechargeable battery with charger.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com