



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Emergency Recovery Trolley

Item Code:

Emergency recovery trolley available in sizes: 210l x cms height 65 to 95 cms. Emergency recovery trolley made of heavy duty erw steel tube, S. S. Swing away railing, S. S. telescopic saline rod, storage tray and oxygen cylinder cage. Emergency recovery trolley are removable x-raytranslucent top and provided with a tray for holding x-ray cassettes underneath the entire length of the stretcher. Height adjustment by hydraulic system, trendelenburg and reverse trendelenburg by pressure cylinder. Fitted with 15 cm dia castors-two with brakes and rubber buffers on all sides.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com