



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Examination Bed

Item Code:

Base frame made of heavy duty square pipe Table top made of urethane foam covered with artificial leather Inclination of Back rest shall be adjustable Detail: Examination Table made of metal frame and metal sheet frame made of square pipe of heavy gauge and top made of sheet well powder coated paint of white color head side adjustable with the manual mechanism as per the requirements of the user provided with rubber shoe. the top of the table fitted with high density foam covered with Rexene cover water proof material overall size of the table 72"x20"x32"

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)