



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Examination Gloves

**Item Code:**

Examination Gloves of free size made of Latex , non sterile packing of 100 pcs powder free, gloves used for external use only for purpose of external examination and prevention, this is not recommended for any surgical procedure.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)