



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Folding Stretcher

Item Code: HF-35

Folding stretcher made of best quality fabric sturdy with folding frame, heavy M.S pipe or Aluminum pipe on the both sides and cross strip provided on the both side of the stretcher to open or to close the stretcher, rubber handle to carry the unit, available in different colors standard size, light weight easily foldable. Fabric quality tetron material with stretcher size 22"x 84" with pole size of 7 fit each side four heavy metal u shape feet provided to hold the stretcher above from the earth. Velcro strap on both side to hold the patient.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com