



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Foot Step Single

Item Code:

Measuring 460 L x 250 W x 230 W mm (single step) single step foot step frame made of M.S pipe round and duly powder coated paint finish, top fixed with anti static material and duly covered with angle on all four side to protect it, the foot step provided with Black rubber shoes anti static type. double foot step also available in the same model.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)