



Gynaecology Examination Table

Item Code:

Gynecological Examination Table. Construction made of stainless steel square or round type tube in two portion of the table one portion duly fitted with wooden top covered with high density of foam and finished with semi leather type material water proof nature, one basin of stainless steel duly fitted with the foot side of this portion of the table that can be easily removable for cleaning , the basin fitted in the stainless steel ring which is also removable two lithotomy legs fitted with foot side of the table and upper portion of the lithotomy legs duly covered with the foam padded and semi leather water proof material for finishing, the lithotomy legs has facility for height adjustment with help of knob, second part of the table which is also constructed with stainless steel and top made of wooden frame with high density of foam and duly covered with semi leather type waterproof material, the part two is movable type and moves on wheels two with lock and two without lock. Dimensions complete table: 180 x 55 x 80 cm.A Delivered with the following accessories: 1 Pair of padded shoulder pads, 1 Stainless bowl, sliding under the seat,1 Pair of padded leg holder.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com