



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Gynecological Examination Table cum Delivery Bed

**Item Code:**

Gynecological examination Table cum Delivery bed made of complete stainless steel construction. The table made of round or Square Tube of stainless steel top of the table fitted with mattress foam padded and duly covered with wash proof semi leather material, head side adjustable with help of the lever provided on the head side, seat part is also adjustable from the lower side, there is movement in the height adjustment of that seat two lithotomy legs provided on the table the upper support for the lithotomy legs duly padded with foam and covered with washable material, velcro provided to lock the legs of the patient on the top of the legs, legs height and movement are adjustable with help of knobs, one basin of stainless steel provided under the U cut of the seat to collect the waste and blood particle during examination and delivery third part of the table are movable towards lower side and can be control with help of handle at different stage of height it is adjustable as per requirements of the user.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)