



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Gynecological Examination Table

Item Code: H.F-25

Gynecological Examination Table made of metal frame of square pipe of heavy gauge head side adjustable with the manual mechanism with the help of handle it allow any desirable position to head side portion as per the requirements of the user. Foot side of the table have mechanism to adjust the height of that portion as per user requirements, there is U cut of the foot side portion and that allow to access the Stainless steel Basin placed on the table. The top of the table fitted with high density foam covered with Rexine cover water proof material overall size of the table 72"x20"x32" the supporting rods with leg holding straps provided with the arrangement of height adjustment as per the requirements, height adjustment of the rods can be placed with the help of the knob. there are head holders on the both side of the table to hold the head of the patient dury the Delivery pains, the head holders are also adjustable as per patient requirements, both the head holder are duly foam padded for comfortable holding of the head.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)