



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Heating Mantle

Item Code: LE-25

Heating Mantle - Heating surface made of glass sleeves reinforced and insulated with glass wool fitted with energy regulator lamp and on/off switch with plug and cord. Work on 220/230V AC.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com