



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Horizontal Autoclave (Round)

Item Code: S.T.R-02

We offer to our clients Horizontal Autoclave (Round) This is a triple jacketed autoclave in Horizontal form. Our range is extensively used in hospitals, private clinics and nursing centers to clean the appliances from the bacterial, fungal and viral infections. Moreover, our range is also extensively used to clean contaminants from the laboratory appliances. This sterilizer disinfects the garbage sacks in buckets that are put in upward order. Internal and external walls of the autoclave are made of high grade stainless steel. Our autoclaves are appreciated for the features like corrosion resistant finish and sturdy construction.

Technical Specifications:

Construction	Sterilizer has inner and outer wall made of 304 grade stainless steel, middle jacket made of stainless steel or brass and gun metal.
Technical Data	Steam Working Pressure 1.26 kgf/cm ² (2.2 kgf/cm ² in case High speed sterilizer)
Working Temperature	121oC
Operating Voltage	240 volts single phase, AC supply, 50 Hz.
Sterilization Period	45 to 50 minutes
Components	Radiant locking system, cold water drainage system
Steam Exhaust	5 to 7 provided with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the pressure inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber made of stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level indicator to see the water inside, heavy duty power plug with socket provided with the unit
Accessories	one spare element, and two gasket , our Regular sizes are 900X500, and 1050X750. Other sizes also available on request.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com