



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT

Hospital Bed Fowler

Item Code: HF-03

Hospital Bed flower that are designed in accordance with High standards. We have designed a strong framework using strong and high quality steel tubes that are mounted on 10 cms diameter on rotating castors for these beds. The top of the perforated CRCA sheets have four sections and are suitably bent with duly machine. Two separate screws are used for the movement of knee rest and backrest. Vertical tubes support is given to the tabular head and foot that are unequal in height. These beds are provided with safe side railings and four saline locations. Furthermore, the finishing is done with epoxy powder coated.

Features:

- High strength
- Sturdy construction
- Size: 198L x 90W x 60H cms

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com