



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Hospital Bed Semi-Fowler

Item Code: HF-06

Hospital Bed semi fowler framework are strong and precise steel tubes having 'H' type legs on rubber stumps. Having two section top of perforated CRCA sheet duly machine press bent with mattress guards, laminated panels fitted with S.S. head & foot bows. Backrest sections maneuvered by crank mechanism from foot end. Provision for four locations for I.V. Rod, safety side railings. Size: 198L x 90W x 60H cms. Finish Epoxy Power Coated.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)