



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Hospital Bed

Item Code:

Hospital Bed main frame are rectangular pipe or Angle Iron, perforated CRCA sheet top tabular head and foot bows of unequal height with vertical tube supports, on rubber stumps. Provision for four saline locations . Size:198L x 90W x 60H cms. Finish: Epoxy Power Coated. available in knockdown model.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com