



HOSPITAL & MEDICAL EQUIPMENT



Hydroculator

Item Code: PHY-13

The three fold fixed upholstered traction table is idle for horizontal cervical and lumbar traction purpose. The lumbar section rolls freely on ball bearing roller for friction free lumbar pull. Cervical traction angle is controlled by vertical adjustment of machine mounting board. Complete with flexion stool. The modified cushioned top traction table is totally folding, which can be dismantled for transportation. Operates in static or intermittent mode. Digital display for hold time and rest time. Feather touch selection for treatment time with digital display. Feather touch selection for treatment time with digital display. Patient safety switch with audio visual alarm. Alarm at the end of treatment comes up and traction force released. Accessories cervical collar, cervical spreader bar, pelvic belt, thoracic belt, traction doubler, main cord.

Technical details:

Power Consumption : 100 watt(approx.) Input voltage : 220 volts a. C. 50 hz. Hold time : 20, 40, 60, 80(second) Release time : 1, 5, 10, 20(second)

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com