



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### ICU Bed Electric

Item Code: HF-01

ICU bed electric A heavy duty, luxury, high finish, electric powered Icu bed is for feature steel frame 1.5mm thickness, scratch free Powder coated paint finish & moulded 4 section bed ventilated mattress base & anti corrosion material motorized backrest motorized knee control motorized hi-low manual support to take over. Battery backup for all movement function by remote Control. trendelenburg & reverse trendelenburg. Foot switch controls hand held control system 125mm castors w/lock special mattress, conductive antibacterial & anti-static. Cornered bumpers / buffer rolls. I. V stands, Side rails on Both side, urine bag holder mattress cover, antibacterial, zipped with liquid proof top & heavy duty.

- Frame work made of rectangular CRCS tubes mounted on 12.5cms dia swiveling castors- two with brakes.
- 18 G perforated CRCA sheet top in four sections.
- Electrically operated Four Separate Actuators for Maneuvering Back rest, knee rest, trendelenburg/reverse trendelenburg & height position by hand remote control system.
- Removable head and foot bows with laminated panels.
- Collapsible safety side railings at both sides.
- S.S. telescopic I.V. rods with four locations.
- Rubber buffers at four corners.
- Pre-treated and Epoxy Powder Coated Finish.
- Size- 218(L) x 98(W) x 55 to 80(H) cms.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : rajasmedical@gmail.com  
www.rajasindia.com