



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Instrument Sterilizer

Item Code: S.T.R-05

A product of choice in the medical field. Seamless shell & lever operated Lid fitted with full proof mechanism control excessive steam escape and restricts condensate within the shell. Incredible synchronized maneuverability of lid, due to statistically perforated tray for flushing & entry of water. Highly adoptable of your growing needs. Configuration design & material as per IS : 5022 Our medical sterilizer, hospital sterilizer available in different sizes: 510 x 200 x 150 mm, 430 x 200 x 150 mm, 300 x 150 x 125 mm. Heating with electric element.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com