



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Instrument Sterilizer

Item Code: S.T.R-05

A product of choice in the medical field. Seamless shell & lever operated Lid fitted with full proof mechanism control excessive steam escape and restricts condensate within the shell. Incredible synchronized maneuverability of lid, due to statistically perforated tray for flushing & entry of water. Highly adoptable of your growing needs. Configuration design & material as per IS : 5022 Our medical sterilizer, hospital sterilizer available in different sizes: 510 x 200 x 150 mm, 430 x 200 x 150 mm, 300 x 150 x 125 mm. Heating with electric element.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)