



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Instruments Cabinet

**Item Code:**

Medical cabinet made of full steel frame with steel angles the cabinet with the provision of the door with glass viewing system from the front, and provided with 3 to 5 shelves inside that is adjustable in height as per the requirements of the user, finished with powder coated paint finish provided with lock and key arrangement, cabinet fitted on the floor feets, shelves of the cabinet made of transparent acrylic sheet or glass as per the demand of the user.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)