



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Instruments Trolley Stainless Steel

Item Code: HF-40

Instruments Trolley size 30x20x20 inches made of tubular pipe complete S.S construction with S.S. top finished units moved on castors Two with breaks and two without breaks provided With two S.S. shelves, handle provided on both sides to move the trolley high edges of shelves from all the four sides to protect the material putting on shelves.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)