



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



IV Stand

Item Code: HF-38

I. V. Stand construction made of steel round type pipe with fully chrome plated on the heavy base unit of five castors to support the whole base, base made of gun metal and black paint epoxy type rod provided totally chrome plated with height adjustable knob on the side to adjust the height of the rod as per the need of the user, two hooks provided on the top of rod to arrange the fluid bottle on the same, provided with bottle holder.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com