



Keratometer

Item Code:

Dual Eye Level Sighting System Dioptre 36 - 52 (in 0.25 steps) Paint Finish : Black Or White The Achromatic Optical Systems to Precise Measurement Of The Corneal Radius Curvature And Accurate Measurement Of Radius Curvature Of Contact Lens.

Keratometer is diagnostic instrument commonly used by ophthalmologists for eye examination, measuring the curvature of the anterior surface of the cornea. Moreover, it provides an accurate measurement of Corneal Radius Curvature and Curvature of Contact lens.

Features:

- Coincidence Focusing System
- 15 x Eyepieces
- One Position Instrument. Measure Both Meridians Without Changing Optical System
- Simple Vertical Adjustable Head To Fit Patient
- Astigmatic Axis 0 to 180 Degree

Standard Accessories:

- Dust Cover
- One Spare Bulb

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com