



HOSPITAL & MEDICAL EQUIPMENT



Lensometer

Item Code:

Lensmeter Corona Cross Target combines an excellent Ocular Type External Reading System with very clear target screen. It is same as Lensmeter Corona Target but with a slight difference with Pinshole Cross + Corona Target. Also the instrument can be moved from 30° to 90° with High Power Magnification for getting sharp and bright corona focus.

Specifications

Type	:	Ocular Type External Reading System
Vertex Power Range	:	± 25D Step : 0.25D Step Up to ± 10D 0.50D Step Beyond ± 10D
Cylindrical Axis Scale	:	0 to 180° (1° step)
Prismatic Power Range	:	5 in the View field (1 step)
Target	:	Pinshole Corona
Eyepiece Focusing	:	+ 3D to -5D
Diameter Of Lens Accommodated	:	30 to 80mm
Tilting Angle	:	30° to 90°
Light Source	:	AC 220V 3W LED

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com