



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Linen Trolley

Item Code: HF29

Linen trolley frame made of round pipe construction unit moves on four wheels of heavy duty type cross pipe bottom frame and four side pipe frame and inner bucket made of thick canvas cloth the canvas bucket upper round duly covered with the cotton rope on his mouth and can be tie once pulled for closing of the mouth of canvas bucket once it filled with gargage or dirty cloths. unit ves on wheels, frame finished with powder coated paint finish, canvas bag dmouly removable.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com