



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Micro Motor Drill

Item Code:

Stepless, smooth speed adjustment from 3,000-35,000 rpm. Precision machining achieved bur eccentricity within 0.02mm. DC coreless micromotor used for quick motor response and smoother rotation over full speed range. The handpiece is compact, lightweight, and well fits to user's hand. User can use it for a long period of time without fatigue. The handpiece is provided with the special dust proof mechanism which prevents entry of cuttings and other foreign substances into the ball bearings, thereby ensuring high durability. Functions and performance desired for laboratory motors are highly emphasized in conception and packaged in a cutely designed control unit. Over-current protector functions immediately and stops the micromotor for safety. Electronically achieved to minimize impact at start to the operator's hand, and this also extends device longevity

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasingia.com