



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Obstetric Delivery Table

Item Code: H.F-28

Obstetric & gynecological divisions Two-section mattress basis with leather or plastic cover suitable for disinfecting Heavy gauge sheet steel top mechanical regulation of back section powder coating paint suitable for wet cleaning Chrome-plated of ratchet supporting cross bar and components. A raising and lowering of a table is actuated through a mechanically (foot paddle) Complete heel stirrups and leg holders with Velcro canvas straps, and leather or plastic covered pad filled with polyurethane or other similar material Support for soles and knees, side protection All steel frame work sturdily braced and mounted on four wheel (two are with lock and two are without lock) Length, approx, mm: 1500 Width, approx, mm 750 Complete Lateral handrail on both sides from stainless steel. Additional Feature I.V. rod for fluids and blood transfusion Basin under top for collection of waste material.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasingia.com](http://www.rajasingia.com)